
Escali

PRODUCT
MANUAL

Ultra Slim Body Composition Scale

Model # BFBW180

www.escali.com

Thank you for choosing the Escali BFBW180 Body Composition Digital Scale. This state-of-the-art, precision measuring instrument is only accurate when used on hard, flat and level surfaces. This scale will not function correctly on carpet. The electronic sensors are very sensitive; avoid dropping or jarring the scale, and please keep it where it will be protected from impact.

Warning! Do not use this product if you have a pacemaker or other internal electronic medical devices.

INITIAL SETUP: CREATING & EDITING USER PROFILES

1. Remove the insulating plastic strip from the battery compartment. Ensure the batteries are properly seated within the +/- terminals.
2. Press the **UNIT** button on the base of the scale to select your desired unit of measure; pounds (lb), kilograms (kg) or stones (st).
3. To set up a new user profile (or edit an existing profile) press **ENTER** and use the **LEFT** or **RIGHT** arrow buttons located on the side of the scale to select a profile number. Press **ENTER** to confirm.
4. "CLrN" will appear on display. Using the arrow buttons on the side, select "N" to continue, or "Y" to reset your profile. Press **ENTER** to confirm.
5. Select your gender using the **LEFT** or **RIGHT** arrow buttons (Activate Athlete Mode by selecting your gender with the accompanied running figure). Press **ENTER** to confirm.
6. Select your height using the **LEFT** or **RIGHT** arrow buttons and press **ENTER** to confirm.
7. Select your age using the **LEFT** or **RIGHT** arrow buttons and press **ENTER** to confirm.
8. When the display shows zero weight, step on the scale with clean and dry bare feet. The scale will measure your weight and then calculate your body composition.

Note: This body scale can save up to 8 unique user profiles. To program additional profiles, follow steps 3 through 8.

SHOULD I ACTIVATE ATHLETE MODE?

Activate Athlete Mode if you exercise at least 10 hours per week and have been doing so consistently for the last 6 months, or longer.

MEASURING BODY WEIGHT & BODY COMPOSITION

1. To use the scale, stand still on the center of the platform with clean and dry bare feet touching each of the metal foot placers. After a few seconds the scale will display your weight.
2. Next, the scale will automatically select your user profile and calculate your body composition while an “o” scrolls across the display.

3. Finally, the scale will display your current body composition measurements. You may now step off the scale. The scale will cycle through your weight, body fat, body water, muscle mass and bone density measurements three times before turning off.

Note: The scale will automatically select your user profile to calculate your body composition by comparing your current body weight measurement to previous body weight measurements from each user profile. If your current body weight changes by more than 7 pounds from your previous measurement, or if your body weight is within 7 pounds of another user profile, then you will need to manually select your user profile number with the corresponding **LEFT** or **RIGHT** arrow buttons before the scale will display your current body composition.

HOW DOES THE SCALE MEASURE BODY COMPOSITION?

This scale uses Bioelectrical Impedance Analysis (BIA) technology to measure body composition. BIA technology sends an unperceivable electronic signal through your body to calculate your body fat, body water and muscle mass percentages. Do not use this scale if you have a pacemaker or other internal electronic medical device.

HOW DO I GET THE MOST ACCURATE RESULTS?

Hydration levels fluctuate throughout the day and are affected by eating, drinking and exercise. As a result, for maximum accuracy and repeatability, it is recommended that this scale is used at the same time of day under the same conditions to establish a personal average. These readings are a close indication, and not an exact value. Please consult your doctor for more information.

MEASUREMENTS DURING WEIGHT LOSS

During weight loss, body fat percentages may not change as expected since the reading is a percentage of body fat compared to body weight, and both could be changing at the same time. To minimize the loss of lean body tissue, be sure to include exercise in your weight loss plan.

BODY FAT PERCENTAGE (5%-60% x 0.1%) — SUGGESTED RANGES

Standard for Men:

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	<13%	<14%	<16%	<17%	<18%
Normal	13-20%	14-21%	16-23%	17-24%	18-25%
Moderately High	20.1-23%	21.1-24%	23.1-26%	24.1-27%	25.1-28%
High	>23%	>24%	>26%	>27%	>28%

Standard for Women:

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	<19%	<20%	<21%	<22%	<23%
Normal	19-28%	20-29%	21-30%	22-31%	23-32%
Moderately High	28.1-31%	29.1-32%	30.1-33%	31.1-34%	32.1-35%
High	>31%	>32%	>33%	>34%	>35%

MUSCLE MASS PERCENTAGE — SUGGESTED RANGES

Men	Approximately 40% of total body weight.
Women	Approximately 30% of total body weight.

BODY WATER PERCENTAGE (43%-73% x 0.1%) — SUGGESTED RANGES

	Actual Body Fat Range	Optimal Body Water Range
Men	4% to 14%	63% to 70%
	14% to 21%	57% to 63%
	22% to 24%	55% to 57%
	25% and over	37% to 55%
Women	4% to 20%	58% to 70%
	21% to 29%	52% to 58%
	30% to 32%	49% to 52%
	33% and over	37% to 49%

BONE MASS PERCENTAGE — SUGGESTED RANGES

The average bone mass percentage for both men and women is between 4% to 5% of total body weight.

PROPER USE & CARE

- Use and stow the scale at room temperature.
- Do not submit the scale to shaking, excessive vibration, or other rough treatment; use on a flat and stable surface.
- The display may be affected by electromagnetic disturbances such as radios, microwaves, or cell phones. If such disturbances occur, remove the source of the disturbance, and restart the scale.
- Disassembling, or tampering with the scale voids the warranty.
- Clean the scale with a damp cloth; do not submerge in liquid; scale is not waterproof.
- Failure to insert batteries in the correct polarity, as indicated in the battery compartment, may shorten the life of the batteries or cause batteries to leak.
- Do not mix Alkaline, Standard (Carbon-Zinc) or Rechargeable (Nickel Cadmium) or (Nickel Metal Hydride) batteries. Do not mix old and new batteries.
- Batteries should be recycled or disposed of as per state and local guidelines. Do not dispose of batteries in fire.

INDICATIONS FOR USE

This device is intended for adults and children (10-17 years old) to measure weight and use bio-electrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, bone mass and muscle mass in a home environment.

TROUBLESHOOTING WARNING MESSAGES

“- - -” — Overload

Solution: Stop using the scale for measurement.

PRODUCT SPECIFICATIONS

Model: BFBW180

Capacity: 400 lb / 180 kg

Increments: 0.2 lb / 0.1 kg

Product Dimensions: 14 x 14 x 1 in / 35 x 35 x 2.4 cm

Power Source: Two CR2032 batteries (included)

ESCALI CORPORATION

3203 Corporate Center Drive, Suite 150

Burnsville, MN 55306

USA

PRODUCT SUPPORT

24/7 Online support: <https://escali.zendesk.com>

Email support: feedback@escali.com

Phone support: 1-800-467-6408

LIFETIME LTD. WARRANTY

This scale is guaranteed to be free of defects in materials and workmanship.

For more details visit: www.escali.com/warranties

Escali's liability is limited to two (2) times the cost of the product.

© Escali Corp., Minneapolis, MN, USA. All rights reserved. Escali, the Escali logo and other Escali marks are owned by Escali Corp. and may be registered.

FCC USER GUIDE INFORMATION

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guar-

antee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna;
- Increase the separation between the equipment and the receiver;
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected;
- Consult the dealer or an experienced radio / TV technician for help.

Radio Transmitters (Part 15)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution

The Escali BFBW180 Ultra Slim Body Composition Bathroom Scale offers you a seamless way to manage your health. Please be aware that this device is designed for adults' self-measuring and self-monitoring body fat and hydration levels. Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening. This device should not be used by anyone who is acutely or chronically ill, suffering from a disease or taking medications that affect your water levels. The accuracy of readings for these patients has not been verified. Specific medical advice should be obtained from a physician.

- Before diet and exercise, you must first find a professional medical guide.
- The Escali BFBW180 is equipped with data transmission function. It may emit electromagnetic energy so as to perform its intended function. Nearby portable and mobile RF communications equipment can affect the performance of the device.
- Portable and mobile RF communications equipment can affect the measuring accuracy of the Escali BFBW180.
- Kindly note that the use of accessories, transducers or cables other than those specified, with internal components, may result in increased EMISSIONS or decreased IMMUNITY of the device.
- Be aware that misuse of electrical equipments can cause electric shock, burns, fire and other hazards. Please do not use the device in a flammable gas environment.
- Warning that the Escali BFBW180 should not be used adjacent to or stacked with other equipment.
- Manufacturer will make available on request circuit diagrams, component parts list, etc.
- WARNING: No modifications of this equipment is allowed.
- During using the scale, the patient will contact with the surface and the electrode of the scale. The surface and the electrode of the scale is type BF Applied part.
- The device doesn't need to be calibrated in five years lifetime.
- Do not touch the output of the batteries when using the analyzer.
- The patient is an intended operator. The patient can perform all the operations in the manual, such as measurement, data transmitting, changing batteries.

Product Name	BFBW180 Ultra Slim Body Composition Scale
Dimension	Scale: 14 x 14 x 1 in / 350 x 350 x 24.5 mm (Approximately)
Net Weight	3.85 lb / 1.75 kg (Approximately)
Display	LCD
Measurement Unit	Pound (lb) / Kilogram (kg) / Stone (st)
Measurement Range	11 lb to 400 lb / 5 kg to 180 kg / 0 st: 11 lb to 28 st: 5 lb
Division	0.2 lb / 0.1 kg / 0 st: 1 lb
Accuracy	50kg: ±0.3kg; 100kg: ±0.4kg; 150kg: ±0.5kg; 180kg: ±0.7kg
Working Environment	Temperature: 41°F to 104°F / 5°C to 40°C; Relative Humidity: ≤90% RH; Atmospheric Pressure: 86 kPa to 106 kPa
Storage Environment	Temperature: -4°F to 140°F / -20°C to 60°C; Relative Humidity: ≤90% RH; Atmospheric Pressure: 50 kPa to 106 kPa
Power Source	6V (Two CR2032-size Batteries)
Turn on Method	Instant On Technology
Auto-OFF	The scale will turn off after about 15 seconds if there is no operation.
Accessories	1. Two CR2032-size batteries; 2. User Manual
Mode of Operation	Continuous Operation
Degree of Protection	Type BF applied part
Software Version	A0
Protection against Ingress of Water	IP21

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "MANUFACTURER"
	Symbol for "TYPE BF APPLIED PARTS"		Symbol for "DIRECT CURRENT"
SN	Symbol for "SERIAL NUMBER"	 	Symbol for "ENVIRONMENT PROTECTION—Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice."
	Symbol for "MANUFACTURE DATE"		

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