# ESCALI BT201 DIGITAL BODY THERMOMETER INSTRUCTIONS


## What is a "normal" temperature?

Although the generally accepted "normal" body temperature is 98.6°F or 37°C, readings can vary from 97.1°F (36.1°C) to 99.0°F (37.2°C) and still may be considered "normal". Variations in temperature can be attributed to activities such as exercise, smoking, eating, drinking and the location of the temperature reading. Also, consider your temperature is often lower in the morning than in the afternoon. While the oral temperature follows the guidelines previously described, a rectal temperature is generally 1.0°F (0.5°C) higher. Conversely, an axillary (under the arm) temperature will usually be 1.0°F (0.5°C) lower.

## Specifications

- Measurement range:** 89.6°F – 109.2°F (32°C – 42.9°C)  
**Low temp. display:** Temperatures <90.0°F (<32.0°C) display L°F (L°C)  
**High temp. display:** Temperatures >109.9°F (42.9°C) display H°F (H°C)  
**Accuracy between 95.9°F – 107.6°F (35.5°C – 42°C):** ±0.2°F (±0.1°C)  
**Display resolution:** 0.1°F (0.1°C)  
**Normal temperature beep sound:** 10 long beep sounds in 10 seconds  
**Fever Alarm sound (>100°F or >37.8°C):** 30 short beep sounds in 10 seconds  
**Auto power off:** 8 minutes (±2 minutes)  
**Battery:** LR41 or AG3  
**Life of battery:** about 1 year or 1200 uses  
**Dimension:** 128 mm x 19 mm x 10 mm  
**Unit Weight:** approx. 9.4 g (including battery)  
**Operating condition:** Temperature 41°F – 104°F (5°C–40°C)  
Humidity 15% – 85% RH  
**Storage condition:** Temperature -4°F – 131°F (-20°C–55°C)  
Humidity 10% – 85% RH

## How to Take your Temperature

1. Press ON/OFF button to activate. The unit will beep and display  for about 2 seconds. (To switch between °F and °C, hold the ON/OFF button for 3 seconds when powering on the thermometer.)
2. After the thermometer powers on, the last temperature taken will be displayed with an "M" in the lower-right corner.
3. When "Lo" displays with °F (°C) flashing, the thermometer is ready for temperature measuring, but if the room temperature is above 90°F (32°C) then room temperature will be displayed instead of "Lo".
4. The beeping will sound when temperature measuring is complete.
5. Press ON/OFF button to power off the unit. (The unit will automatically power off approximately eight minutes after completing the measurement.)

### Taking an Oral Temperature Measurement

Place the probe tip under the tongue, breath through the nose and keep the lips closed. Approximate measuring time is 1 minute.

### **Taking an Axillary Temperature (Under the Arm)**

Wipe the underarm with a dry towel and place the probe tip under the arm so the tip is touching the skin with the thermometer perpendicular to the body. Position the arm across the chest so the probe tip is well covered by the arm.

Approximate measuring time: 2~3 minutes

Note: Axillary temperatures are generally 1.0°F (0.5°C) lower than oral readings.

### **Taking a Rectal Temperature**

Apply a water-soluble lubricant to the probe cover. Gently insert the probe into the rectum (a maximum 1/2 inch.)

Approximate measuring time: 1 minute

Note: Rectal temperatures are generally 1.0°F (0.5°C) higher than oral readings.

## **Replacing the Battery**

The battery in your digital thermometer needs to be replaced when “ $\uparrow$ ” appears in the lower-right corner of the display. Replace the battery with a 1.5V LR41 silver oxide or AG3 alkaline battery (or equivalent.)

### **To replace the battery:**


1. Pull the battery cover away from the unit.

2. Pull out the battery holder about 18mm.

3. Use suitable tool to remove the old battery and add a new battery.

## **Cleaning Instructions**

Before and after each use, clean the thermometer with a soft cloth and isopropyl alcohol diluted with water, or wash with cold soapy water. Do not soak for an extended period of time or sterilize by boiling method, gas or steam autoclaves.

## **General Care**

1. The thermometer should only be used under the supervision of an adult.
2. Do not walk, run, eat, drink or talk during temperature measurement.
3. Store the unit in the protective case when not in use.
4. Do not store the unit where it will be exposed to direct sunlight, dust or humidity. Avoid extreme temperatures.
5. Dropping or subjecting your thermometer to strong shocks should be avoided.
6. Do not attempt to disassemble the unit, except to replace the battery.

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