
Complete Health Body Composition Scale

Model # BF180-2


www.escali.com

Thank you for choosing the Escali BF180-2 Complete Health Body Composition Scale. This state-of-the-art, precision measuring instrument is only accurate when used on hard, flat and level surfaces. This scale will not function correctly on carpet. The electronic sensors are very sensitive; avoid dropping or jarring the scale, and please keep it where it will be protected from impact.

Warning! Do not use this product if you have a pacemaker or other internal electronic medical devices.

DISPLAY DIAGRAM

	Body Fat Analysis Results		Male
	Total Body Water Analysis Results		Female
	Muscle Mass Analysis Results		Male Athlete
	Body Bone Analysis Results		Female Athlete
CAL	Calorie	%	Percentage
BMI	Body Mass Index	AGE	Age
ft	Foot	st lb	Stones : Pounds
cm	Centimeter	lb	Pounds
++	Obese	kg	Kilograms
+	Over Fat	▢	Under Fat
INITIAL	Healthy	SETUP:	CREATING & EDITING USER PROFILES

1. Remove the insulating plastic strip from the battery compartment. Ensure the batteries are properly seated within the +/- terminals.
2. Press the **UNIT** button on the base of the scale to select your desired unit of measure; pounds (lb), kilograms (kg) or stones (st : lb).
3. To set up a new user profile (or edit an existing profile) press **SET** and use the **UP**

or **DOWN** arrows to select a profile number. Press **SET** to confirm.

4. "CLrN" will appear on display. Using the **UP** or **DOWN** arrows, select "N" to continue, or "Y" to reset your profile. Press **SET** to confirm.
5. Select your gender using the **UP** or **DOWN** arrows. (Activate Athlete Mode by selecting your gender with the accompanied running figure.) Press **SET** to confirm.
6. Select your height using the **UP** or **DOWN** arrows and press **SET** to confirm.
7. Select your age using the **UP** or **DOWN** arrows and press **SET** to confirm.
8. When the display shows zero weight, step on the scale with clean and dry bare feet. The scale will measure your weight and then calculate your body composition.

Note: This body scale can save up to 8 unique user profiles. To program additional profiles, follow steps 3 through 8.

SHOULD I ACTIVATE ATHLETE MODE?


Activate Athlete Mode if you exercise at least 10 hours per week and have been doing so consistently for the last 6 months, or longer.

MEASURING BODY WEIGHT & BODY COMPOSITION

1. To use the scale, stand still on the center of the platform with clean and dry bare feet touching each of the metal foot placers. After a few seconds the scale will display your weight.
2. Next, the scale will automatically select your user profile and calculate your body composition while a "o" scrolls across the display.

3. Finally, the scale will display your current body composition measurements. You may now step off the scale. The scale will cycle through your body fat, body water, muscle mass and bone density measurements three times before

itself off.

Note: The scale will automatically select your user profile to calculate your body composition by comparing your current body weight measurement to previous body weight measurements from each user profile. If your current body weight changes by more

than 4.5 pounds from your previous measurement, or if your body weight is within 4.5 pounds of another user profile, then you will need to manually select your user profile number. For example, if profile 1 and 2 are within 4.5 pounds the scale will display "P1P2" on the LCD. You can select the correct profile number by pressing the **UP** or **DOWN** buttons under your profile number. The scale will proceed to display body composition measurements.

HOW DOES THE SCALE MEASURE BODY COMPOSITION?

This scale uses Bioelectrical Impedance Analysis (BIA) technology to measure body composition. BIA technology sends an unperceivable electronic signal through your body to calculate your body fat, body water and muscle mass percentages. Do not use this scale if you have a pacemaker or other internal electronic medical device.

HOW DO I GET THE MOST ACCURATE RESULTS?

Hydration levels fluctuate throughout the day and are affected by eating, drinking and exercise. As a result, for maximum accuracy and repeatability, it is recommended that this scale is used at the same time of day under the same conditions to establish a personal average. These readings are a close indication, and not an exact value. Please consult your doctor for more information.

MEASUREMENTS DURING WEIGHT LOSS

During weight loss, body fat percentages may not change as expected since the reading is a percentage of body fat compared to body weight, and both could be changing at the same time. To minimize the loss of lean body tissue, be sure to include exercise in your weight loss plan.

BODY FAT PERCENTAGE — SUGGESTED RANGES

Men—

Age	Under fat	Healthy	Over fat	Obese
≤ 29	< 13%	13-20%	20.1-23%	> 23%
30-39	< 14%	14-21%	21.1-24%	> 24%
40-49	< 16%	16-23%	23.1-26%	> 26%
50-59	< 17%	17-24%	24.1-27%	> 27%
60 ≥	< 18%	18-25%	25.1-28%	> 28%

Women—

Age	Under fat	Healthy	Over fat	Obese
≤ 29	< 19%	19-28%	28.1-31%	> 31%
30-39	< 20%	20-29%	29.1-32%	> 32%
40-49	< 21%	21-30%	30.1-33%	> 33%
50-59	< 22%	22-31%	31.1-34%	> 34%
60 ≥	< 23%	23-32%	32.1-35%	> 35%

MUSCLE MASS PERCENT— SUGGESTED

AGE

RANGES

Age	6-12	13-20	21-40	41-60	61-80
Men	> 39%	> 37%	> 36%	> 35%	> 33%
Women	> 35%	> 33%	> 32%	> 31%	> 30%

TOTAL BODY WATER (TBW) PERCENTAGE — SUGGESTED RANGES

	Body Fat % Range	Optimal TBW % Range
	Men	4 to 14%
15 to 21%		63 to 57%
22 to 24%		57 to 55%
25% and over		55 to 37%
Women	4 to 20%	70 to 58%
	21 to 29%	58 to 52%
	30 to 32%	52 to 49%
	33% and over	49 to 37%

BONE MASS (LB/KG) — SUGGESTED RANGES

Body Weight	Men			Women		
	< 132.3 lb < 60 kg	132.3 - 165.3 lb 60 - 75 kg	> 165.3 lb > 75 kg	< 99.1 lb < 45 kg	99.1 - 132.3 lb 45 - 60 kg	> 132.3 lb > 60 kg
Est. Bone Mass	5.5 lb 2.5 kg	6.5 lb 2.9 kg	7 lb 3.2 kg	4 lb 1.8 kg	5 lb 2.2 kg	5.5 lb 2.5 kg

PROPER USE & CARE

- Use and stow the scale at room temperature.
- Do not submit the scale to shaking, excessive vibration, or other rough treatment; use on a flat and stable surface.
- The display may be affected by electromagnetic disturbances such as radios, microwaves, or cell phones. If such disturbances occur, remove the source of the disturbance, and restart the scale.
- Disassembling, or tampering with the scale voids the warranty.
- Clean the scale with a damp cloth; do not submerge in liquid; scale is not water-proof.
- Failure to insert batteries in the correct polarity, as indicated in the battery compartment, may shorten the life of the batteries or cause batteries to leak.
- Do not mix Alkaline, Standard (Carbon-Zinc) or Rechargeable (Nickel Cadmium) or

(Nickel Metal Hydride) batteries. Do not mix old and new batteries.

- Batteries should be recycled or disposed of as per state and local guidelines. Do not dispose of batteries in fire.

TROUBLESHOOTING WARNING MESSAGES

“- - -” — Overload

Solution: Stop using the scale for measurement.

“Lo”— Low Battery

Solution: Replace both batteries as measurements may lose accuracy.

“P1P2” — Unable to automatically select a user profile

Solution: Use the UP or DOWN arrow buttons to select your user profile.

PRODUCT SPECIFICATIONS

Model: BF180-2

Capacity: 400 lb / 180 kg

Increments: 0.2 lb / 0.1 kg

Product Dimensions: 11.8 x 11.8 x 0.94 in (30 x 30 x 2.4 cm)

Power Source: Two CR2032 batteries (included)

ESCALI CORPORATION

3203 Corporate Center Drive, Suite 150

Burnsville, MN 55306, USA

PRODUCT SUPPORT

24/7 Online support: <https://escali.zendesk.com>

Email support: feedback@escali.com

Phone support: 1-800-467-6408

LIFETIME LTD. WARRANTY

This scale is guaranteed to be free of defects in materials and workmanship. For more details visit: www.escali.com/warranties

Escali's liability is limited to two (2) times the cost of the product.

© Escali Corp., Minneapolis, MN, USA. All rights reserved. Escali, the Escali logo and other Escali marks are owned by Escali Corp. and may be registered.

FCC USER GUIDE INFORMATION

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates,

uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna;
- Increase the separation between the equipment and the receiver;
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected;
- Consult the dealer or an experienced radio / TV technician for help.

Radio Transmitters (Part 15)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution

The Escali BF180-2 Complete Health Body Composition Scale offers you a seamless way to manage your health. Please be aware that this device is designed for adults' self-measuring and self-monitoring body fat and hydration levels. Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening. This device should not be used by anyone who is acutely or chronically ill, suffering from a disease or taking medications that affect your water levels. The accuracy of readings for these patients has not been verified. Specific medical advice should be obtained from a physician.

- Before diet and exercise, you must first find a professional medical guide.
- The Escali BF180-2 is equipped with data transmission function. It may emit electromagnetic energy so as to perform its intended function. Nearby portable and mobile RF communications equipment can affect the performance of the device.
- Portable and mobile RF communications equipment can affect the measuring accuracy of the Escali BF180-2.
- Kindly note that the use of accessories, transducers or cables other than those specified, with internal components, may result in increased EMISSIONS or decreased IMMUNITY of the device.
- Be aware that misuse of electrical equipments can cause electric shock, burns, fire and other hazards. Please do not use the device in a flammable gas environment.
- Warning that the Escali BF180-2 should not be used adjacent to or stacked with other equipment.
- Manufacturer will make available on request circuit diagrams, component parts list, etc.
- WARNING: No modifications of this equipment is allowed.
- During using the scale, the patient will contact with the surface and the electrode of the scale. The surface and the electrode of the scale is type BF Applied part.
- The device doesn't need to be calibrated in five years lifetime.
- Do not touch the output of the batteries when using the analyzer.
- The patient is an intended operator. The patient can perform all the operations in the manual, such as measurement, data transmitting, changing batteries.

EMC GUIDANCE

Table 1 —

Product Name	BF180-2 Complete Health Body Composition Scale
Dimension	Scale: 11.8 x 11.8 x 0.94 in / 300 x 300 x 24 mm (Approximately)
Net Weight	3.41 lb / 1.55 kg (Approximately)
Display	LCD
Measurement Unit	Pound (lb) / Kilogram (kg) / Stone (st)
Measurement Range	11 lb to 397 lb / 5 kg to 180 kg / 0 st: 11 lb to 28 st: 5 lb
Division	0.2 lb / 0.1 kg / 0 st: 1 lb
Accuracy	50kg: ±0.3kg; 100kg: ±0.4kg; 150kg: ±0.5kg; 180kg: ±0.7kg
Working Environment	Temperature: 50°F to 104°F / 10°C to 40°C; Relative Humidity: 15% RH to 90% RH; Atmospheric Pressure: 80 kPa to 105 kPa
Storage Environment	Temperature: -4°F to 131°F / -20°C to 55°C; Relative Humidity: ≤90% RH
Power Source	6V (Two CR2032-size Batteries)
Turn on Method	Instant On Technology
Auto-OFF	The scale will turn off after about 15 seconds if there is no operation.
Accessories	1. Two CR2032-size batteries; 2. User Manual
Mode of Operation	Continuous Operation
Degree of Protection	Type BF applied part
Software Version	A0
Protection against Ingress of Water	IP21

	Symbol for “THE OPERATION GUIDE MUST BE READ”		
	Symbol for “TYPE BF APPLIED PARTS”		Symbol for “DIRECT CURRENT”
SN	Symbol for “SERIAL NUMBER”		Symbol for “ENVIRONMENT PROTECTION—Waste electrical products should not be disposed of with household
	Symbol for “MANUFACTURE DATE”		Please recycle where facilities exist. Check with your local authority or retail-

Guidance and Manufacturer's Declaration — Electromagnetic Emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not application
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not application

Table 2 —

Guidance and Manufacturer's Declaration — Electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Not application	Not application
Surge IEC 61000-4-5	Not application	Not application
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not application	Not application
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	Not application	Not application
Radiated RF IEC61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
NOTE U _T is the a.c. mains voltage prior to application of the test level.		

Table 3 —

Guidance and Manufacturer's Declaration - Electromagnetic Immunity							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
	870						
	930						
	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
5240							
5785							

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